

Anti-Rat IgG3 antibody (STJ16102550)

GENERAL INFORMATION

| | |
|---------------|---------------------|
| Product Type | Primary antibodies |
| Applications | FACS/IF/IHC-F/IHC-P |
| Host / Source | Rat |
| Reactivity | Rat |

PRODUCT PROPERTIES

| | |
|---------------------|----------------------------------|
| Clonality | Polyclonal |
| Concentration | 0.5 mg/mL |
| Conjugation | Unconjugated |
| Purification | Affinity Purified |
| Formulation | TRIS, PBS, pH7 and antibiotics |
| Isotype | IgG3 |
| Storage Instruction | Store at -20°C for up to 1 year. |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|