

Anti-CD25/IL2RA antibody [143-13] (STJ16102559)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/FC/FUNC/WB/IF/IHC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	ELISA (solid phase: 0, 1-100 µg/ml; tracer: 0, 001-100 µg/ml for 20 min at RT). Flow cytometry (1-2 µg/million cells for 30 min, at 4°C). Functional studies (0, 02-2, 0 µg/ml without azide). Immunoblotting (1-2 µg/ml). Immunofluorescence (1-2 µg/ml).
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG1k
Storage Instruction	Store for up to 1 year at 2-8°C upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---