

Anti-CD36 antibody [185-1G2] (STJ16102561)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IF/IHC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	Flow cytometry (1-2 µg/million cells in 0, 1 ml for 30 min, at 4°C). Immunofluorescence (1-2 µg/ml). Immunohistology (1-2 µg/ml for 30 min at RT; no antigen retrieval procedure is known to date for staining of formalin-fixed tissues).
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG2bk
Storage Instruction	Store for up to 1 year at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	948
Gene Symbol	CD36
UniProt ID	CD36_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---