

Anti-Estrogen receptor antibody [AER-314] (STJ16102569)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC
Host / Source	Mouse
Reactivity	Cow/Human/Mouse/Rabbit/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	Immunohistology (formalin-fixed: 2-4 µg/ml for 30 min at RT; staining of formalin-fixed tissues requires boiling tissues sections in 10mM citrate buffer, pH 6.0, for 10-20 min followed by cooling at RT for 20 minutes).
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG1 Lambda
Storage Instruction	Store for up to 1 year at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	2852
Gene Symbol	GPER1
UniProt ID	GPER1_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---