

Anti-CD 5 antibody [BB8] (STJ16102576)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | FC/IF |
| Host / Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Clonality | Monoclonal |
| Concentration | Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL. |
| Conjugation | Unconjugated |
| Purification | Affinity purified from tissue culture. |
| Dilution Range | Flow cytometry (0, 5-1, 0 µg/million cells in 0, 1 ml). Immunofluorescence (0, 5-1, 0 µg/ml). |
| Formulation | PBS with 0.02% Sodium Azide. |
| Isotype | IgG1k |
| Storage Instruction | Store for up to 1 year at 2-8°C upon receipt. |

TARGET INFORMATION

| | |
|-------------|---------------------------|
| Gene ID | 921 |
| Gene Symbol | CD5 |
| UniProt ID | CD5_HUMAN |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|