

Anti-Keratin 19 antibody [E6] (STJ16102581)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IF/WB/IHC
Host / Source	Mouse
Reactivity	Rat/Rabbit/X.laevis

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	Flow Cytometry (1-2 µg/million cells in 0.1 ml, fix cells in 4% PFA for 10 min, at 4°C, permeabilize with 0, 2% saponin or digitonin for 15 min, at 4°C). Immunofluorescence (1-2 µg/ml). Immunoblotting (1-2 µg/ml). Immunohistology (1-2 µg/ml for 30 mi
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG1k
Storage Instruction	Store for up to 1 year at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	360626
Gene Symbol	Krt19
UniProt ID	K1C19_RAT

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---