

Anti-NuMA antibody [EBS-C-002] (STJ16102582)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/IHC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	ELISA (solid phase: 0, 1-100 µg/ml; tracer: 0, 001-100 µg/ml for 30 min at RT). Immunoblotting (1-2 µg/ml). Immunohistology (1-2 µg/ml for 30-60 minutes at RT; staining of formalin-fixed tissues requires boiling tissue sections in 10mM Tris with 1mM
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgGmK
Storage Instruction	Store for up to 1 year at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	4926
Gene Symbol	NUMA1
UniProt ID	NUMA1_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---