

Anti-CD2 antibody [EBS-CD-002] (STJ16102584)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/WB/IF/IHC
Host / Source	Mouse
Reactivity	Cat/Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	Flow cytometry (0.5-1.0 µg/million cells in 0.1 ml). Immunoblotting (1.0-2.0 µg/ml). Immunofluorescence (0.5-1.0 µg/ml). Immunohistology (1.0-2.0 µg/ml for 30 min at RT; an appropriate antigen retrieval method for staining of formalin-fixed tissues h
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG2ak
Storage Instruction	Store for up to 1 year at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	914
Gene Symbol	CD2
UniProt ID	CD2_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---