

Anti-CD19 antibody [EBS-CD-015] (STJ16102585)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IF/IHC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	Flow Cytometry (0, 5-2, 0 Mu g/million cells in 0, 1ml for 30 min, at 4°C). Immunofluorescence (0, 5-2, 0 Mu g/ml). Immunohistology (1-2 Mu g/ml for 30 min at RT; no suitable antigen retrieval procedure is known to date for formalin-fixed tissues).
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG1k
Storage Instruction	Store for up to 1 year at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	12478
Gene Symbol	Cd19
UniProt ID	CD19_MOUSE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---