

Anti-Keratin 18 antibody [EBS-IF-006] (STJ16102597)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	FC/WB/IF/IHC
	Host / Source	Mouse
	Reactivity	Cat/Dog/Human/Pig
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
	Conjugation	Unconjugated
	Purification	Affinity purified from tissue culture.
	Dilution Range	Flow Cytometry (1-2 µg/million cells in 0.1 ml, fix cells in 4% PFA for 10 min, at 4°C, permeabilize with 0, 2% saponin or digitonin for 15 min, at 4°C). Immunoblotting (1-2 µg/ml). Immunofluorescence (1-2 µg/ml). Immunohistology (formalin-fixed: 2-4
	Formulation	PBS with 0.02% Sodium Azide.
	Isotype	IgG1k
	Storage Instruction	Store for up to 1 year at 2-8°C upon receipt.
TARGET INFORMATION		
	Gene ID	3875
	Gene Symbol	KRT18
	UniProt ID	K1C18_HUMAN
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.