

Anti-TNF-alpha antibody [EBS-O-174] (STJ16102602)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/FC/WB/IF/IHC
Host / Source	Mouse
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	ELISA (solid phase: 0, 1-100 µg/ml; tracer: 0, 001-100 µg/ml for 30 min at RT). Flow cytometry (0, 5-1, 0 ug/million cells in 0.1ml). Immunoblotting (1-2 µg/ml). Immunofluorescence (1-2 Mu g/ml). Inhibits tumor growth (order Ab without azide). Immuno
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG1 Kappa
Storage Instruction	Store for up to 1 year at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	7124
Gene Symbol	TNF
UniProt ID	TNFA_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---