

Anti-Neurofilaments, myomesin-NF-H antibody [NE14] (STJ16102624)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB/IHC
Host / Source	Mouse
Reactivity	Cat/Chicken/Cow/Gerbil/Guinea Pig/Human/Mouse/Pig/Rabbit/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	Immunoblotting (0, 5-1, 0 µg/ml). Immunohistology (1-2 µg/ml for 30 min at RT; staining of formalin-fixed tissues requires boiling tissue sections in 10mM citrate buffer, pH 6.0, for 10-20 min followed by cooling at RT for 20 minutes).
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG1k
Storage Instruction	Store for up to 1 year at 2-8°C upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---