

## Anti-CD 1b antibody [RIV12] (STJ16102627)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/FC/IF/IHC
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	ELISA (solid phase: not known; tracer: 0.001-100 µg/ml for 30 min at RT). Flow cytometry (0.5-1.0 µg/million cells in 0.1 ml). Immunofluorescence (0.5-1.0 µg/ml). Immunohistology (1.0-2.0 µg/ml for 30 min at RT; an appropriate antigen retrieval metho
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG1k
Storage Instruction	Store for up to 1 year at 2-8°C upon receipt.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------