

Anti-Keratin 17 antibody [E3] (STJ16102660)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IHC-F/IHC-P/IF/WB
Host / Source	Mouse
Reactivity	Cat/Cow/Dog/Goat/Human/Pig/Rabbit/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 Mu g/mL
Conjugation	Unconjugated
Formulation	PBS with 0.02% sodium azide
Isotype	IgG2bk
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---