

**Anti-Keratin 4, 5, 6, 8, 10, 13, 18 antibody [C-11] (STJ16102663)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC/IHC-F/IHC-P/IF/WB
Host / Source	Mouse
Reactivity	Cat/Chicken/Cow/Dog/Ferret/Goat/Guinea Pig/Human/Kangaroo/Marmoset/Monkey/Mouse/Rabbit/Rat/Xenopus Laevis

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	100 Mu g/mL
Conjugation	Unconjugated
Formulation	PBS with 0.02% sodium azide
Isotype	IgG1k
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---