

Anti-Vimentin antibody [LN-6] (STJ16102664)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IHC-F/IHC-P/IF
Host / Source	Mouse
Reactivity	Cat/Cow/Human/Mouse/Pig/Rabbit/Rat/Sheep

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 Mu g/mL
Conjugation	Unconjugated
Formulation	PBS with 0.02% sodium azide
Isotype	IgMk
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	7431
Gene Symbol	VIM
UniProt ID	VIME_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---