

Anti-Campylobacter jejuni antibody [EBS-I-002] (STJ16102682)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IF/IHC-F
Host / Source	Mouse
Reactivity	C. Jejuni/Type1/Type2

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 Mu g/mL
Conjugation	Unconjugated
Formulation	PBS with 0.02% sodium azide
Isotype	IgMk
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---