

Anti-Clostridium difficile Toxin A antibody [EBS-I-100] (STJ16102683)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IHC-F/IF
Host / Source	Mouse
Reactivity	C. Difficile

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 Mu g/mL
Conjugation	Unconjugated
Formulation	PBS with 0.02% sodium azide
Isotype	IgG3k
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------