

Anti-MRP14/S100A9/calgranulin B antibody [47-8D3] (STJ16102707)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC/IF/IHC-F/IHC-P/WB
Host / Source	Mouse
Reactivity	Baboon/Cat/Cow/Dog/Goat/Guinea Pig/Human/Horse/Monkey/Pig/Rabbit/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 Mu g/mL
Conjugation	Unconjugated
Formulation	PBS with 0.02% sodium azide
Isotype	IgG1k
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---