

## Anti-CD19 antibody [CB19] (STJ16102720)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IHC-F/IF
Host / Source	Mouse
Reactivity	Human/Mouse

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 Mu g/mL
Conjugation	Unconjugated
Formulation	PBS with 0.02% sodium azide
Isotype	IgG1k
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### TARGET INFORMATION

Gene ID	<a href="#">930</a>
Gene Symbol	<a href="#">CD19</a>
UniProt ID	<a href="#">CD19_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---