

Anti-CD4 antibody [EDU-2] (STJ16102732)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IHC-F/IF
Host / Source	Mouse
Reactivity	Human/Chimpanzee/Orangutan/Stump-Tailed Macaque

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 Mu g/mL
Conjugation	Unconjugated
Formulation	PBS with 0.02% sodium azide
Isotype	IgG2ak
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	920
Gene Symbol	CD4
UniProt ID	CD4_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---