

Anti-Laminin antibody (STJ16102775)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/IHC-P/WB
Host / Source	Rabbit
Reactivity	Bovine/Human/Mouse/Pig/Rabbit/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Dilution Range	IHCF IHC-P 1:50-1:100 WB
Formulation	Contains 0.09% sodium azide
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---