

Anti-VMAT2 antibody (STJ16102790)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F
Host / Source	Rabbit
Reactivity	Dog/Guinea Pig/Hamster/Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Dilution Range	IHCF 1:800-1:1, 200
Formulation	Lyophilized; reconstitute in 100 µl dist. water
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	6571
Gene Symbol	SLC18A2
UniProt ID	VMAT2_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---