

Anti-Amylin antibody (STJ16102791)

GENERAL INFORMATION

| | |
|---------------|------------------------------|
| Product Type | Primary antibodies |
| Applications | IHC-F/IHC-P |
| Host / Source | Rabbit |
| Reactivity | Cat/Guinea Pig/Human/Pig/Rat |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Polyclonal |
| Conjugation | Unconjugated |
| Dilution Range | IHCF IHC-P 1:800-1:1, 500 |
| Formulation | Lyophilized; reconstitute in 100 μ l dist. water |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt. |

TARGET INFORMATION

| | |
|-------------|----------------------------|
| Gene ID | 3375 |
| Gene Symbol | IAPP |
| UniProt ID | IAPP_HUMAN |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|