

Goat Anti-Mouse IgG/A/M H&L antibody {FITC} (STJ16102802)

GENERAL INFORMATION

Product Type	Secondary antibodies
Applications	FC
Host / Source	Goat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	1mg/mL
Conjugation	FITC
Purification	Purified
Formulation	Phosphate buffered saline<0.1% Sodium Azide (NaN ₃)
Storage Instruction	Store at-20°C for up to 1 year. Minimise freeze-thaw cycles.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---