

Anti-COLLAGEN IV antibody (STJ16102806)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IB/IF/IHC-F/ELISA
Host / Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.4 mg/mL
Conjugation	Unconjugated
Purification	Purified
Formulation	Borate buffered salineNone present
Storage Instruction	Store at-20°C for up to 1 year. Minimise freeze-thaw cycles.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---