

Goat Anti-Human IgG/A/M antibody {FITC} (STJ16102811)

GENERAL INFORMATION

Product Type	Secondary antibodies
Applications	ELISA/FC
Host / Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	2.0 mg/mL
Conjugation	FITC
Purification	lyophilized
Formulation	Phosphate buffered saline0.09% Sodium Azide (NaN ₃)
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---