

Anti-Serotonin (5-HT) antibody (STJ16102814)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/IHC-P/IF
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Dilution Range	IHCF IHC-P IF 1:200-1:400
Formulation	Lyophilized; reconstitute in 100 μ l dist. water
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---