

## Anti-Helospectin antibody (STJ16102818)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IHC-F              |
| Host / Source | Rabbit             |
| Reactivity    | Cat/Dog/Mouse/Rat  |

### PRODUCT PROPERTIES

|                     |                                                 |
|---------------------|-------------------------------------------------|
| Clonality           | Polyclonal                                      |
| Conjugation         | Unconjugated                                    |
| Dilution Range      | IHCF 1:400-1:800                                |
| Formulation         | Lyophilized; reconstitute in 100 µl dist. water |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt.  |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|