

Anti-Insulin/Proinsulin antibody (STJ16102819)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/IHC-P/IF
Host / Source	Guinea pig
Reactivity	Hamster/Human/Pig/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Dilution Range	IHCF IHC-P 1:3, 000-1:5, 000 IF
Formulation	Lyophilized; reconstitute in 100 µl dist. water
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---