

## Anti-Pancreastatin antibody (STJ16102820)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IHC-F/IHC-P        |
| Host / Source | Rabbit             |
| Reactivity    | Rat                |

### PRODUCT PROPERTIES

|                     |                                                 |
|---------------------|-------------------------------------------------|
| Clonality           | Polyclonal                                      |
| Conjugation         | Unconjugated                                    |
| Dilution Range      | IHCF<br>IHC-P 1:200-1:1, 000                    |
| Formulation         | Lyophilized; reconstitute in 100 µl dist. water |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt.  |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|