

## Anti-Histamine antibody (STJ16102821)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/IHC-P/IF
Host / Source	Rabbit
Reactivity	Human/Rat

### PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Dilution Range	IHCF IHC-P 1:1, 500-1:2, 000 IF
Formulation	Lyophilized; reconstitute in 100 $\mu$ l dist. water
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---