

## Anti-Glucagon antibody (STJ16102825)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/IHC-P
Host / Source	Guinea pig
Reactivity	Cat/Human

### PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Dilution Range	IHC-F IHC-P 1:1, 000-1:2, 000
Formulation	Contains 0.09% sodium azide
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---