

Anti-VIP antibody (STJ16102826)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/IHC-P
Host / Source	Guinea pig
Reactivity	Cat/Guinea Pig/Human/Pig/Rabbit/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Dilution Range	IHCF IHC-P 1:1, 000-1:2, 000
Formulation	Lyophilized; reconstitute in 100 µl dist. water
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	7432
Gene Symbol	VIP
UniProt ID	VIP_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---