

Anti-PYY antibody (STJ16102830)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/IHC-P/IF
Host / Source	Guinea pig
Reactivity	Cat/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Dilution Range	IHCF IHC-P IF 1:1, 500-1:2, 500
Formulation	Lyophilized; reconstitute in 100 μ l dist. water
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	5697
Gene Symbol	PYY
UniProt ID	PYY_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---