

Anti-Adrenomedullin antibody (STJ16102832)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F
Host / Source	Guinea pig
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Dilution Range	IHCF 1:400-1:800
Formulation	Lyophilized; reconstitute in 100 µl dist. water
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---