

Anti-NIS, Na⁺/I Symporter antibody (STJ16102840)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/IF
Host / Source	Rabbit
Reactivity	Human/Pig/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Dilution Range	IHCF IF 1:1, 200-1:2, 000
Formulation	Lyophilized; reconstitute in 100 µl dist. water
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---