

Anti-Human IgM-Fc antibody (STJ16102842)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IP
Host / Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Formulation	lyophilized, add 1 ml sterile distilled water
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------