

Anti-SOX10 antibody [SQab1734] (STJ16102864)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/ICC/IF/IHC-P/IP/WB
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Purification with Protein A.
Formulation	PBS, 0.01% Sodium azide, 40% Glycerol and 0.05% BSA.
Storage Instruction	Store at -20°C for up to 1 year. Minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6663
Gene Symbol	SOX10
UniProt ID	SOX10_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---