

## Anti-ERK2 antibody [SQab1744] (STJ16102874)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FC/IHC-P/IP/WB     |
| Host / Source | Rabbit             |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |  |
|---------------------|--|
| Clonality           | Monoclonal   |
| Conjugation         | Unconjugated   |
| Purification        | Purification with Protein A.                                 |
| Formulation         | PBS, 0.01% Sodium azide, 40% Glycerol and 0.05% BSA.         |
| Storage Instruction | Store at-20°C for up to 1 year. Minimise freeze-thaw cycles. |

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|