

## Anti-CDK4 antibody [SQab1864] (STJ16102878)

### GENERAL INFORMATION

|               |                       |
|---------------|-----------------------|
| Product Type  | Primary antibodies    |
| Applications  | FC/ICC/IF/IHC-P/IP/WB |
| Host / Source | Rabbit                |
| Reactivity    | Human                 |

### PRODUCT PROPERTIES

|                     |                                                               |
|---------------------|---------------------------------------------------------------|
| Clonality           | Monoclonal                                                    |
| Conjugation         | Unconjugated                                                  |
| Purification        | Purification with Protein A.                                  |
| Formulation         | PBS, 0.01% Sodium azide, 40% Glycerol and 0.05% BSA.          |
| Storage Instruction | Store at -20°C for up to 1 year. Minimise freeze-thaw cycles. |

### TARGET INFORMATION

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">1019</a>       |
| Gene Symbol | <a href="#">CDK4</a>       |
| UniProt ID  | <a href="#">CDK4_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|