

## Anti-HDAC1 antibody [SQab1867] (STJ16102880)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/ICC/IF/IHC-P/WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Purification with Protein A.
Formulation	PBS, 0.01% Sodium azide, 40% Glycerol and 0.05% BSA.
Storage Instruction	Store at -20°C for up to 1 year. Minimise freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">3065</a>
Gene Symbol	<a href="#">HDAC1</a>
UniProt ID	<a href="#">HDAC1_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---