

**Anti-Cytokeratin 19 antibody [SQab1872] (STJ16102885)****GENERAL INFORMATION**

|               |                       |
|---------------|-----------------------|
| Product Type  | Primary antibodies    |
| Applications  | FC/ICC/IF/IHC-P/IP/WB |
| Host / Source | Rabbit                |
| Reactivity    | Human                 |

**PRODUCT PROPERTIES**

|                     |                                                               |
|---------------------|---------------------------------------------------------------|
| Clonality           | Monoclonal                                                    |
| Conjugation         | Unconjugated                                                  |
| Purification        | Purification with Protein A.                                  |
| Formulation         | PBS, 0.01% Sodium azide, 40% Glycerol and 0.05% BSA.          |
| Storage Instruction | Store at -20°C for up to 1 year. Minimise freeze-thaw cycles. |

**ADDITIONAL INFORMATION**

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|