

## Anti-PAX2 antibody [SQab1896] (STJ16102902)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P
Host / Source	Rabbit
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Purification with Protein A.
Formulation	PBS, 0.01% Sodium azide, 40% Glycerol and 0.05% BSA.
Storage Instruction	Store at-20°C for up to 1 year. Minimise freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">5076</a>
Gene Symbol	<a href="#">PAX2</a>
UniProt ID	<a href="#">PAX2_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---