

## Anti-Cytokeratin 5 antibody [SQab18106] (STJ16102915)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IHC-P/WB           |
| Host / Source | Rabbit             |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |                                                               |
|---------------------|---------------------------------------------------------------|
| Clonality           | Monoclonal                                                    |
| Conjugation         | Unconjugated                                                  |
| Purification        | Purification with Protein A.                                  |
| Formulation         | PBS, 0.01% Sodium azide, 40% Glycerol and 0.05% BSA.          |
| Storage Instruction | Store at -20°C for up to 1 year. Minimise freeze-thaw cycles. |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|