

Anti-p21 antibody [SQab18118] (STJ16102923)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/ICC/IF/IHC-P
Host / Source	Rabbit
Reactivity	Other

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Purification with Protein A.
Formulation	PBS, 0.01% Sodium azide, 40% Glycerol and 0.05% BSA.
Storage Instruction	Store at -20°C for up to 1 year. Minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1026
Gene Symbol	CDKN1A
UniProt ID	CDN1A_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---