

Anti-Cytokeratin 8 + 18 antibody [SQab20229] (STJ16102983)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	IHC-P
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Purification with Protein A.
Formulation	PBS, 0.01% Sodium azide, 40% Glycerol and 0.05% BSA.
Storage Instruction	Store at -20°C for up to 1 year. Minimise freeze-thaw cycles.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---