

Anti-c-Myc antibody [SQab22253] (STJ16103002)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | IHC-P |
| Host / Source | Rabbit |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Monoclonal |
| Conjugation | Unconjugated |
| Purification | Purification with Protein A. |
| Formulation | PBS, 0.01% Sodium azide, 40% Glycerol and 0.05%BSA. |
| Storage Instruction | Store at-20°C for up to 1 year. Minimise freeze-thaw cycles. |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|