

**Anti-PHOX2B antibody [SQab22261] (STJ16103009)****GENERAL INFORMATION**

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IHC-P              |
| Host / Source | Rabbit             |
| Reactivity    | Human              |

**PRODUCT PROPERTIES**

|                     |                                                              |
|---------------------|--------------------------------------------------------------|
| Clonality           | Monoclonal                                                   |
| Conjugation         | Unconjugated                                                 |
| Purification        | Purification with Protein A.                                 |
| Formulation         | PBS, 0.01% Sodium azide, 40% Glycerol and 0.05%BSA.          |
| Storage Instruction | Store at-20°C for up to 1 year. Minimise freeze-thaw cycles. |

**TARGET INFORMATION**

|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">8929</a>        |
| Gene Symbol | <a href="#">PHOX2B</a>      |
| UniProt ID  | <a href="#">PHX2B_HUMAN</a> |

**ADDITIONAL INFORMATION**

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|