

Anti-ROR gamma antibody [SQab30332] (STJ16103043)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	IHC-P
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Purification with Protein A.
Formulation	PBS, 0.01% Sodium azide, 40% Glycerol and 0.05%BSA.
Storage Instruction	Store at-20°C for up to 1 year. Minimise freeze-thaw cycles.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---